

Hiden HPR-20 QIC TMS Transient MS for Fast Event Gas Analysis



Introduction

The Hiden HPR-20 QIC TMS is configured for continuous analysis of gases and vapours at pressures near atmosphere.

Operating to 200°C, the QIC (quartz inert capillary) flexible 1m capillary inlet provides fast response times of less than 150 ms.

The HPR-20 QIC TMS system has a mass range of 200 AMU (300, 510, 1000 AMU options) and a detection capability from 100% to less than 5 ppb.

The TMS system is optimised for fast response studies and can respond to changes in gas composition with a 5 decade response time in <200 ms.

TMS – for Fast Event Studies

Features of the system to optimise the response for fast event studies include:

- Pulse Ion Counting detection with 7 decade continuous log scale
- 1m capillary gives <150 ms response time
- Open ion source increases pumping speed on ion source



Ultra-fast response



HPR-20 QIC TMS: > 5 decades response in < 0.2 s

Ultra-fast response



The HPR-20 QIC TMS tracks the changes in concentration of fast gas pulses with incredible response and accuracy over more than 5 decades. This data shows the measurement of 5 pulses in 5 seconds.

HPR-20 TMS Vacuum Schematic



Key

G1 Penning gauge
VR1 QIC Inlet bypass control valve
P1 60 I/s turbo drag pump
P2 Backing and bypass Scroll pump
MS UHV Housing (Mass spectrometer chamber)



Backing and bypass Scroll Pump

HPR-20 QIC TMS Analyser:

Hiden HAL 3F PIC Triple Filter Mass Spectrometer



Triple Filter Mass Spectrometer

Why have a triple filter?

Two main advantages:

- 1. Strict control over the quadrupole entrance and exit fields provides enhanced sensitivity for high mass transmission and increased abundance sensitivity
- 2. Enhanced long-term stability. The bulk of the deselected ions from the quadrupole ioniser deposit harmlessly on the RF-only pre-filter stage, minimising contamination on the mass selective primary filter.



QIC Inlet Technology



 Quartz and Platinum Wetted Surfaces
 →
 No memory effects

 Heated Capillary
 →
 No condensation effects

 Flow Matched
 →
 Optimum response / recovery

 Minimal Internal Volume
 →
 PPB detection

 Interchangeable Sampling Capillaries
 →
 Analysis from 10 Torr to 2 Bar

Typical Mass Spectrum of Air



Note: Different species can have the same mass e.g. CO, N_2 m/e 28

Soft Ionisation

Unique to Hiden gas analysis systems, soft ionisation allows users to selectively ionise different gases by setting the ionisation energy for a particular mass.

This powerful technique can simplify the analysis of otherwise complex cracking patterns from multi-component gas/vapour mixtures.

The ionisation energy can be altered from 4 to 150 eV, in 0.1 eV increments. Standard operation is at 70 eV.



/ Water

🔚 HA... 😰 🔍 🖾

/ Oxygen

/ Carbon Dioxide

14789 - Serie

MASsoft 7 Professional control software



HPR-20 QIC TMS - Applications

- Fast event studies
- Reaction kinetics
- Selective Catalytic Reduction (SCR)
- Steady State Isotopic Transient Kinetic Analysis (SSITKA)
- High speed switching analysis
- Operando studies
- SpaciMS



Applications: Catalysis Research

A kinetic study of the effect of H_2 on the Selective Catalytic Reduction of NO_x with octane using isotopically labelled ¹⁵NO, using an HPR-20 QIC TMS system.



Figure 1. Mass spectrometer response for ${}^{15}N_2$, ${}^{15}NO$, H_2 and Kr when switching 0.72% H_2 in and out of a SCR feed stream over the catalyst at 300 °C.

Ref: J. P. Breen, R. Burch, C. Hardacre, C. J. Hill and C. Rioche (2007). *J. Catal.*, 2007, 246, 1, p1-9.

Applications: Breath Analysis

Breath by Breath Analysis of Expired Isoprene during Exercise



HPR-20 QIC TMS data showing breath by breath isoprene levels during an exercise test.

SAMSUNG





Hiden HPR-20 Users





UNIVERSITY OF CAMBRIDGE

Imperial College London





Eidgenössische Technische Hochschule Zürich Swiss Federal Institute of Technology Zurich

NASA **Dow Chemical** Exxon-Mobil Imperial College MIT University of British Columbia University of Queensland BASF Seoul National University Suzuki University of Cambridge **Beijing Institute of Technology** Samsung ETH Zürich KAUST **Durham University** Siemens Shell



Summary

- Specifically designed for fast event studies less than 150 ms response time
- Bench-top triple filter quadrupole mass spectrometer gas analysis system
- Real-time, multi-species analysis 5 PPB to 100%
- Soft ionisation for reduced spectral fragmentation and simplified data interpretation





Hiden Analytical Ltd. 420 Europa Boulevard Warrington, WA5 7UN, England

www.HidenAnalytical.com

info@hiden.co.uk

Tel: +44 (0)1925 445 225

