

Hiden HPR-20 QIC R&D

for Advanced Research

Quadrupole Mass Spectrometers for Advanced Science




Introduction

The Hiden HPR-20 QIC R&D is configured for continuous analysis of gases and vapours at pressures near atmosphere.

Operating to 200°C, the QIC (quartz inert capillary) flexible 2m capillary inlet provides fast response times of less than 300 ms.

The HPR-20 QIC R&D system has a mass range of 200 AMU (300, 510, 1000 AMU options) and a detection capability from 100% to less than 5 ppb.

HPR-20 Vacuum Schematic




Key

G1 Penning gauge VR1 QIC Inlet bypass control valve P1 60 l/s turbo drag pump P2 Backing and bypass Scroll pump MS UHV Housing (Mass spectrometer chamber)



Backing and bypass Scroll Pump

HPR-20 QIC R&D Analyser: Hiden HAL 3F Series Triple Filter Mass Spectrometer




Triple Filter Mass Spectrometer

Why have a triple filter?

Two main advantages:

- 1. Strict control over the quadrupole entrance and exit fields provides enhanced sensitivity for high mass transmission and increased abundance sensitivity
- 2. Enhanced long-term stability. The bulk of the deselected ions from the quadrupole ioniser deposit harmlessly on the RF-only pre-filter stage, minimising contamination on the mass selective primary filter.

QIC Inlet Technology



Quartz and Platinum Wetted Surfaces No memory effects


Heated Capillary No condensation effects

Flow Matched Optimum response / recovery

Minimal Internal Volume PPB detection


Interchangeable Sampling Capillaries ——— Analysis from 10 Torr to 2 Bar

Fast Response to Permanent Gases / Vapours



Data shows the response of a HPR-20 system to gas and vapour during switching between a dry He stream and a wet H_2 and Ar flow. For clarity, only the H_2 and H_2 O data is shown in the graph.

Typical Mass Spectrum of Air



Note: Different species can have the same mass e.g. CO, N_2 m/e 28

Soft Ionisation

Unique to Hiden gas analysis systems, soft ionisation allows users to selectively ionise different gases by setting the ionisation energy for a particular mass.

This powerful technique can simplify the analysis of otherwise complex cracking patterns from multi-component gas/vapour mixtures.

The ionisation energy can be altered from 4 to 150 eV, in 0.1 eV increments. Standard operation is at 70 eV.





Figure 1 A: m/z vs Electron energy-H₂O/Air


Figure 1 B: m/z vs Electron energy-NH₃/ H₂O/Air mix


Gas Sampling Options

Pressure:

Inlet options are available for sampling both above and below atmospheric pressure. High pressure inlets for sampling at up to 30 bar and special capillaries for sampling down to 1 mbar.








Multi-stream selectors:

2, 8, 16, 20, 40 and 80 way options


Temperature:

- Heated capillary extensions
- High temperature capillary inlets
- Hot-zone adaptors
- Heated multi-stream inlets

MASsoft 7 Professional control software



QGA Professional Software for Quantitative Gas Analysis



An application specific software package for quantitative gas and vapour analysis providing real time continuous analysis of up to 32 species with concentrations measured in the range 0.1PPM to 100%.

- Automatic subtraction of spectral overlaps
- Automated calibration routines
- Mass spectral library with intelligent scan feature
- Multi-stream support


Applications

- Catalysis
- Reaction Kinetics
- TPD/TPR/TPO
- Thermal Analysis Mass Spectrometry
- Gas Purity Analysis
- Process Characterisation
- Fermentation Off Gas Analysis
- Environmental Gas Analysis
- Combustion Studies
- CVD/MOCVD



Applications: Catalysis Research

- Catalyst characterisation
- Kinetic and thermodynamic measurements
- TPD, TPO, TPR, TP-Reaction
- On-line continuous product analysis
- Total Surface Area / Metal Surface Area
- Mechanisms of Surface Reactions
- Heats of Adsorption and Co-adsorption
- Operando Studies




 NH_3 -TPD on embedded Ru@LSZ and impregnated Ru/LSZ catalysts.

Ref: Lorenzut et al. (2011) *Hydrogen production* through alcohol steam reforming on Cu/ZnO-based catalysts, Applied Catalysis B, **101** (3&4), 397-408.

Applications: Environmental Gas Analysis







Breakthrough curve of CO_2 (15% CO_2 , 85% balance N_2) on mesoporous alumina.

Ref: Yang et al. (2010) CO_2 adsorption over ionexchanged zeolite beta with alkali and alkaline earth metal ions, Mesoporous Materials **135** (1-3), 90-94.

Quadrupole Mass Spectrometers for Advanced Science





Hiden HPR-20 Users





NASA **Dow Chemical** Exxon-Mobil Imperial College MIT

BASF

Seoul National University

University of Cambridge

Samsung ETH Zürich **KAUST**

Durham University Siemens



University of British Columbia University of Queensland

Suzuki

Beijing Institute of Technology

Shell







Massachusetts Institute of Technology













Imperial College London







Eidgenössische Technische Hochschule Zürich Swiss Federal Institute of Technology Zurich

Summary

- Bench-top triple filter quadrupole mass spectrometer gas analysis system
- Real-time, multi-species analysis 5 PPB to 100%
- Fast response to permanent gases and vapours – less than 300 ms response time
- Soft ionisation for reduced spectral fragmentation and simplified data interpretation



Quadrupole Mass Spectrometers for Advanced Science



Hiden Analytical Ltd. 420 Europa Boulevard Warrington, WA5 7UN, England

www.HidenAnalytical.com

info@hiden.co.uk

Tel: +44 (0)1925 445 225

